


FMT

Medical Sensors, Cables and Accessories



ECG CABLES FOR PATIENT MONITORING





Yoke Type can be determined from [Cable Yokes](#) section and you get detailed information about yoke types.

Compatible Brand/Model can be seen at [Compatibility Table](#) . You can find the part number of the cable for desired brand & model with number of leadwires.

Patient End connection types can be seen at the end of the [Compatibility Table](#).


Color Code of the cable can be selected between IEC or AHA color coding options at [Color Codes and Placement](#) section.

Replacement Leadwires are available for FMT or other ECG cable yokes. Please see the [Leadwires](#) Section for compatible leadwires.

Patient End connection types; Snap (S), Grabber (G), Veterinary Clip (V).



• Snap



• Grabber



• Veterinary Clip

ECG Cables for Patient Monitoring

Cable Yokes



• E200-5xxx

E200 yoke for one piece cable with fixed leads.
Yoke cable 3.0 m, lead wire 0.70 m



• E200-4xxx

E200 yoke for one piece cable with fixed leads
Yoke cable 3.0 m, lead wire 0.70 m



• E200-3xxx

E200 yoke for one piece cable with fixed leads
Yoke cable 3.0 m, lead wire 0.70 m



• E201-5xxx

E201 yoke for safety DIN 1.6 mm leadwire.
Complete cable comes with E21/5X or E22/5X leadwires as default. Yoke cable 3.0 m, lead wire 0.70 m



• E201-3xxx

E201 yoke for safety DIN 1.6 mm leadwire.
Complete cable comes with E31/5X or E22/3X leadwires as default. Yoke cable 3.0 m, lead wire 0.70 m



• E202-5xxx
• E202N-5xxx

E202 yoke for with dual pole Siemens style grouped leadwire
Datascope cable yokes
Complete cable comes with E31/5X or E31/5XN and E32/5X or E32/5XN leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m



• E202-3xxx

E202 yoke for dual pole Siemens style grouped leadwire
Complete cable comes with E31/3X or E32/3X leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m



• E204-3xxx
• E204-5xxx

E204 yoke for dual pole GE style grouped leadwire.
Complete cable comes with E41/5XT or E42/5XT and E41/3XT or E42/3XT leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m





• E205-5xxx

E205 yoke for dual pole HP style one piece grouped leadwire.
Complete cable comes with E51/5X or E52/5X leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m.



• E205-3xxx

E205 yoke for dual pole HP style one piece grouped leadwire.
Complete cable comes with E51/3X or E52/3X leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m.



• E206-5xxx

E206 yoke for dual pole Spacelabs style leadwire.
Complete cable comes with E61/5X or E62/5X leadwires as default. Yoke cable 3.0 m, lead wire 0.70 m



• E206-3xxx

E206 yoke for dual pole Spacelabs style leadwire.
Complete cable comes with E61/3X or E62/3X leadwires as default. Yoke cable 3.0 m, lead wire 0.70 m



• E207-3xxx

• E207-6xxx

E207 yoke for Nihon Kohden style one piece grouped leadwire.
Complete cable comes with E71/3X and E71/6X or E72/3X and E726X leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m



• E208-5xxx

E208 yoke for Philips one piece grouped leadwire.
Complete cable comes with E81/5X or E82/5X leadwires as default. Yoke cable 3.0 m, lead wire 0.95 m



• E208-3xxx

E208 yoke for Philips one piece grouped leadwire.
Complete cable comes with E81/3X or E82/3X leadwires as default.



• E209-6xxx

E209 yoke for one piece cable with fixed leads.
yoke cable 1.44 m, lead wire 0.62 m

ECG Cables for Patient Monitoring

Cable Yokes



• E209-5xxx

E209 yoke for one piece cable with fixed leads
Yoke cable 1.44 m, lead wire 0.62 m



• E209-4xxx

E209 yoke for one piece cable with fixed leads
Yoke cable 1.44 m, lead wire 0.62 m





• E209-3xxx

E209 yoke for one piece cable with fixed leads
Yoke cable 1.44 m, lead wire 0.62 m



Color Codes and Placement



ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Adecon



E20X-3001/YZ
E20X-5001/YZ

● DK-8000/D/L/M/S.

● Connector No: 001



E20X-3001/YZ
E20X-5001/YZ

● iM 7/8, PM 900/900S, BM 9000/9000S.

● Connector No: 001



E20X-3007/YZ
E20X-5007/YZ

● AnyView A2E/3/5/6/8, EMS, Q3/5/7.

● Connector No: 036



E20X-3001/YZ
E20X-5001/YZ

● M8000/9000.

● Connector No: 001



E20X-3077/YZ
E20X-5077/YZ

● M7000/8500/9500.

● Connector No: 022



E20X-3051/YZ

● BM3.

● Connector No: 051



E20X-3050/YZ
E20X-5050/YZ

● BM5, BM7.

● Connector No: 049



E20X-3053/YZ

● BPM103.

● Connector No: 140



E20X-3001/YZ
E20X-5001/YZ

● C80, STAR8000A/B/C/D/H.

● Connector No: 001

Biocare

Biolight

Bionet

Biosys

Comen

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Comen



E20X-3031/YZ ● C100, C70, C50, C60, C90, C20, C30, NC8.
E20X-5031/YZ

● Connector No: 248



E20X-3094/YZ ● BP88NET/NEXT/SNEXT, BP306, BP508.
E20X-5094/YZ

● Connector No: 076



E20X-3093/YZ ● BP88SINEXT, BPS510.
E20X-5093/YZ

● Connector No: 066



E20X-4058/YZ ● Corpuls 08/16.

● Connector No: 058



E209-4059/YZ ● Corpuls 1, Corpuls 3

● Connector No: 266



E209-6059/YZ ● Corpuls 1, Corpuls 3
(Complementary ECG Diagnostic Cable)

● Connector No: 267

Criticare Systems (CSI)



E20X-3001/YZ ● 507, 507E, Scholar EL, Scholar II 507EP, Scholar III 507EL,
E20X-5001/YZ 508, 1100, 8100E nGenuity, 8100EL nGenuity, 8100H nCompass, 8500 Poet IQ2.

● Connector No: 001



E20X-3002/YZ ● 2000/A, 2001/A, 2100, 2200I, 3000/A, M/D3 w/Leadselect,
E20X-5002/YZ M/D3A, M/D4, Passport 2/LT, Passport EL/XG/5L, Baloon Pump System 90/90T, 97, 98, Spectrum, Trio, Passport II.

● Connector No: 001

Datascope



E20X-5028/YZ ● DX2022, DX2515.

● Connector No: 039

Dixtal



E20X-3030/YZ ● Sirecust 400 / 500 / 600 / 700 / 900 / 200 / 341 / 610 / 620 / 630 /
E20X-5030/YZ 720 / 722 / 730 / 732 / 960 / 961 / 1280 / 1281 / 1801, Q960.

● Connector No: 030

Dräger Siemens

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Dräger Siemens



● Connector No: 083



E20X-3083/YZ

● Dräger VM1/3, RM1, PM8010/8014, Babyguard 8000, Cato, Babylog 8000/8000 plus, PB8800, Cicero.



● Connector No: 174



E201-3098/GZ

● Neomed Pod Multiparameter Cable.

Emtel



● Connector No: 022



E20X-3048/YZ

● FX2000.

E20X-5048/YZ



● Connector No: 163



E20X-3047/YZ

● FX3000.

E20X-5047/YZ

Fukuda Denshi



● Connector No: 067



E20X-3014/YZ

● DS5100E, DS5300W, DS5400, DS7001.

E20X-5014/YZ



● Connector No: 043



E20X-3043/YZ

● DS7000 / 7100 / 7200 / 7300.

E20X-5043/YZ

GE Datex



● Connector No: 010



E20X-3010/YZ

● AS/3, AS/3 Compact, CS/3, S/5 (with Datex Module), S/5 Light, Cardiocap, Cardiocap II, Cardiocap 5, Light, Satlite Plus.

E20X-5010/YZ

GE



● Connector No: 060



E20X-3060/YZ

● B20 / 40, S/5, Eagle 1000 / 3000 / 5000, Solar 8000 / 9500, Dash 1000 / 2500 / 4000 / 5000, TRAM 100 / 200 / 300 / 400 / 500, E-PRESTN, E-PRETN, PDM module, MAC-LAB.

E204-3060/YZ(*)

E20X-5060/YZ

E204-5060/YZ(**)

(*) Provided with E4X-3XT leadwire set. (**) E4X-5X and E4X-5XT leadwire sets are available for E204 cable yokes. Please indicate your selection in brackets after the part number. EXAMPLE: E204-5060/SI (E42/SI)

Goldway Philips



● Connector No: 024



E20X-3024/YZ

E205-3024/YZ

E208-3024/YZ

E20X-5024/YZ

E205-5024/YZ

E208-5024/YZ

● G60/70/80.



● Connector No: 001



E20X-3001/YZ

● UT4000 A/F/A Pro/F Pro, UT6000/F, G30/40.

E20X-5001/YZ

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Hellige



E20X-3069/YZ
E20X-5069/YZ

● Eagle 1000, Cardioserv, SMV and SMV series, Vicom-SM, Servomed4.

● Connector No: 069



E20X-3091/YZ

● The LifePulse 110 R.

● Connector No: 022



E20X-3040/YZ
E20X-5040/YZ

● Smartsigns Compact 750/1000.

● Connector No: 157



E20X-5006/YZ

● Innocare-S, Innocare-T12.

● Connector No: 148



E20X-3001/YZ
E20X-5001/YZ

● MDE Escort M8/M10/M12, Millenia, Omni Trak, Escort/II, Prism.

● Connector No: 001

Invivo



E20X-3022/YZ
E20X-5022/YZ

● Arrow KAAT II Plus Intra Aortic Balloon Pump, 7135, 7135B, 7141, Kolormon, Micromon, Minimon, Supermon, Vitalmon.

● Connector No: 022

Kontron



E20X-3042/YZ

● LM.

● Connector No: 042

Medtronic Physio Control



E200-3057/YZ
E200-5057/YZ

● LifePak 12/15/20/20e.

● Connector No: 243

Mek



E20X-3089/YZ
E20X-5089/YZ

● MP500.

● Connector No: 089

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Mek



E20X-3079/YZ
E20X-5079/YZ

- MP700/800(Pamo II).

● Connector No: 080



E20X-3029/YZ
E20X-5029/YZ

- MP1000NT, MP1300.

● Connector No: 029



E200-3073/YZ

- Melacard.

● Connector No: 150



E20X-3037/YZ
E20X-5037/YZ

- Horizon 900 / 1100 / 1110 / 2000 / 2200 / 5200, 9000XL, MR1200, 1270, 1280, 1300, 1380.

● Connector No: 037



E20X-3075/YZ
E20X-5075/YZ

- Envoy, Menntor X7, Vitalogik 4000 / 4500 / 5000 / 5500 / 6000.

● Connector No: 075

Menna



E200-4015/YZ

- Defi N/B, Defimonitor ECO1, DM1/3/10/30.

● Connector No: 015



E200-4065/YZ

- Defimonitor ECO1, DM1/3.

● Connector No: 065



E200-3033/YZ

- Heartsave 6/6S, Defimonitor XD.

● Connector No: 033



E200-4016/YZ

- Defimonitor XD, XD1, XD3, XD300.

● Connector No: 169

Metrax Primedic

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Metrax Primedic



E200-3066/YZ
E200-5066/YZ

● Metrax Primedic Defimonitor EVO, Mediana D500.

● Connector No: 262



E200-5001/YZ

● Life Point (Monophasic).

● Connector No: 001



E200-5071/YZ

● Life Point.

● Connector No: 133



E200-4027/YZ

● Life Point (Biphasic).

● Connector No: 098

Metsis

Mindray



E20X-3001/YZ
E20X-5001/YZ

● PM7000/8000/9000, MEC 1000/1200.

● Connector No: 001



E20X-3084/YZ
E20X-5084/YZ

● PM5000/6000, iPM 8/10/12/9800, iMEC 8/10/12, BeneView T1/T5/T6/T8/T9.

● Connector No: 084

Monitor



E20X-3011/YZ
E20X-5011/YZ

● Mitar.

● Connector No: 157

Nec



E20X-3045/YZ
E20X-5045/YZ

● 3M12.

● Connector No: 245

Nihon Kohden



E200-3009/YZ

● Cardiolife TEC-5521 / 5531 / 7100 / 7200 / 7621 / 7631 / 7721 / 7731.

● Connector No: 009

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Nihon Kohden



● Connector No: 090

E20X-3090/YZ

E207-3090/YZ

E207-6090/YZ

● LifeScope i BSM-2301 / 2303 / 2304, LifeScope L BSM-2153 / 2353 / 2351, LifeScope P BSM-4101 / 4102 / 4103 / 4104 / 4111 / 4112 / 4113 / 4114, LifeScope A BSM-5105 / 5135, LifeScope TR BSM-6301 / 6501 / 6701, LifeScope J BSM-9101, LifeScope M 9510, LifeScope VS BSM-3500 / 3700, LifeScope PT BSM-1700, PVM 2701, Vismo PVM 2703, CardioLife TEC-8321 / 8322 / 8332 / 8342 / 8352.



● Connector No: 204

E20X-3086/YZ

E207-3086/YZ

● LifeScope N OPV 1500K.



● Connector No: 008

E20X-3008/YZ

E20X-5008/YZ

● LifeScope 6 OEC-6102 / 6105 / 6120.



● Connector No: 001

E20X-3001/YZ

E20X-5001/YZ

● LifeScope 8 OEC-8108, LifeScope 9 BSM-8301 / 8302 / 8500 / 8502, LifeScope 6 OEC-6102 / 6105 / 7102 / 8108, CardioLife TEC-5200 / 7100 / 7200 / 7300, TEC- 7511 / 7521 / 7531 / 7541.

Pace-Tech



● Connector No: 001

E20X-3001/YZ

E20X-5001/YZ

● Minipack 300 / 3000 / 3100, Vitalmax 800 / 2100 / 2200 / 4000 / 4100.

Palco



● Connector No: 004

E20X-5044/YZ

● 500, 500VP.

Petaş



● Connector No: 140

E20X-3003/YZ

E200-4003/YZ

E20X-5003/YZ

● PM100, KMA275/375/675/800.



● Connector No: 004

E20X-4004/YZ

E20X-5004/YZ

● KMA 165/265/365.

Philips HP



● Connector No: 020

E20X-3020/YZ

E205-3020/YZ

E20X-5020/YZ

E205-5020/YZ

● 43100A / 43110A/M / 43120A/MC, 43200A, M3500B, 7700AC / 77020A/ 77010AC, 77025A, 77030A, 77035A, 77040A / 78208A / 78210A / 7330A / 78331A / 78332A / 78333A / 78341A, 78342A / 78345A / 78341A / 78351T / 78352A/C / 78353A/B / 78354A / 78532A/B / 78354A/B/C / 78619A / 78620A / 78660A/B / 78661A/B / 78670A / 78671A / M1722A/B / M2400A / M2406A / Donor 2000/2500 / M2425A / M2475B.

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

 <p>Connector No: 024</p>		<p>E20X-3024/YZ E205-3024/YZ E208-3024/YZ E20X-5024/YZ E205-5024/YZ E208-5024/YZ</p>	<p>Philips HP</p> <ul style="list-style-type: none"> 78352C, 78354C, 78834C, 862474(C3), M1001A, M1001B, M1002A, M1165A, M1166A, M1167A, M1175A, M1177A, M1205A, M1264A, M1275A, M1280A, M1723A/B, M2475B, M2601A, M3000A, M30001A, M3046A, M4M3, M4735A.
------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

 <p>Connector No: 022</p>		<p>E200-3035/YZ E200-5035/YZ</p>	<p>Progetti</p> <ul style="list-style-type: none"> Rescue LIFE defibrillator.
------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	--------------------------------------	-----------------------------------------------------------------------------------------------------

 <p>Connector No: 096</p>		<p>E200-4097/YZ</p>	<ul style="list-style-type: none"> Rescue 230.
------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	---------------------	---------------------------------------------------------------

 <p>Connector No: 022</p>		<p>E20X-3092/YZ E20X-5092/YZ</p>	<p>RGB</p> <ul style="list-style-type: none"> Omicrom FT Surveyor, Omicrom Vision.
------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	--------------------------------------	----------------------------------------------------------------------------------------------------------

 <p>Connector No: 085</p>		<p>E20X-3085/YZ E20X-5085/YZ</p>	<p>Saadat</p> <ul style="list-style-type: none"> Novin S630/S740/S1600/S1800, Zagros.
--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------	--------------------------------------	-------------------------------------------------------------------------------------------------------------

 <p>Connector No: 022</p>		<p>E20X-3081/YZ E20X-4081/YZ E20X-5081/YZ</p>	<p>Schiller Bruker - Odam</p> <ul style="list-style-type: none"> Argus LCM / LCM Plus / TM7, CM8,
--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------	-------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------

 <p>Connector No: 022</p>		<p>E200-3082/YZ E200-4082/YZ</p>	<ul style="list-style-type: none"> Defiguard 2002/3002/4000/5000/6002, SM784/785, Physiogard 2000/2005/3002/TM 910.
--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------	--------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

 <p>Connector No: 038</p>		<p>E20X-3038/YZ E206-3038/YZ E20X-5038/YZ E206-5038/YZ</p>	<p>Spacelabs</p> <ul style="list-style-type: none"> Ultraview 90309 / 90496 / 90367 / 90369 / 90385 / 90387 / 91496 / 91415 / SL2200 / 2400 / 2600 / 2700, qube, élance, élance elite, Xprezzon, mCare 300 91220.
--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------	------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

 <p>Connector No: 001</p>		<p>E20X-3001/YZ E20X-5001/YZ</p>	<ul style="list-style-type: none"> 90300 / 90303 / 90305 / 90308 / 90309 / 90369 / 90385 / 90408 / 90418 / 90425 / 90428 / 90432 / 90470 / 90700 Series, TEK-208 / 408 / 413 / 414 / 511 / 514 / 521.
--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------	--------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

ECG Cables for Patient Monitoring

Compatible Brands & Models

▼ Connector

▼ Part Number

▼ Brand / Model

Spengler



E20X-3070/YZ
E20X-5070/YZ

● Monica 2000.

● Connector No: 153



E20X-3010/YZ
E20X-5010/YZ

● Diascope 1/2, Diascope NT, Athena, Cardio-Aid, Traveller, Neoscope.

S&W (Artema)

● Connector No: 010



E20X-5096/YZ

● UM300.

Utas

● Connector No: 096



E20X-3001/YZ
E20X-5001/YZ

● UM300-10 / 12 / 15 / 500.

● Connector No: 001



E20X-3001/YZ
E20X-5001/YZ

● Visio.

Üzümcü

● Connector No: 001



E20X-3007/YZ
E20X-5007/YZ

● Visio Mod.

● Connector No: 036



E20X-3012/YZ

● Micropaq/Propaq LT. Vital Sign Monitor.

Wellch Allyn

● Connector No: 048



E20X-3002/YZ
E20X-5002/YZ

● Propaq 100/200 Series, Encore, CS; Ultra Smart.

● Connector No: 001



E200-3002/YZ
E200-5002/YZ

● D900, PD1200/1400/1600.

Zoll


● Connector No: 001

ECG Cables for Patient Monitoring

Leadwires

▼ Yoke Side

▼ Patient Side



- E21-5A ● E21-5I
- E21-3A ● E21-3I

DIN type leadwires, 0.70 m (for E201 cable yokes)



- E22-5I ● E22-5A
- E22-3I ● E22-3A



- E31-3I ● E31-3A

SIEMENS/DRÄGER leadwire, 0.95 m (for E202-3XXX cable yokes)



- E32-3A ● E32-3I



- E32-3I-OR ● E32-3A-OR

SIEMENS/DRÄGER leadwire with ESU protection (Internal ECG ESU BLOCK), 0.95 m (for E202-3XXX cable yokes)

Only Grabber.



- E31-5I ● E31-5A

SIEMENS/DRÄGER leadwire, 0.95 m (for E202-5XXX cable yokes)



- E32-5I ● E32-5A



- E32-5I-OR ● E32-5A-OR

SIEMENS/DRÄGER leadwire with ESU protection (Internal ECG ESU BLOCK), 0.95 m (for E202-5XXX cable yokes)

Only Grabber.



- E31-3IM ● E31-3AM

3 leadwire setting for Dräger cable yokes with 2 dummy leads.

SIEMENS/DRÄGER leadwire, 0.95 m (for E202-5XXX cable yokes)



- E32-3IM ● E32-3AM



- E32-3IM-OR ● E32-3AM-OR

3 leadwire setting for Dräger cable yokes with 2 dummy leads.

SIEMENS/DRÄGER leadwire with ESU protection (Internal ECG ESU BLOCK), 0.95 m (for E202-5XXX cable yokes)

Only Grabber.



- E31-5IN ● E31-5AN

5 leadwire setting for Datascope cable yokes.

SIEMENS/DRÄGER leadwire, 0.95 m (for E202-5XXX cable yokes)



- E32-5IN ● E32-5AN

● IEC color code ● AHA color code

ECG Cables for Patient Monitoring

Leadwires

▼ Yoke Side

▼ Patient Side



● E31-3IN ● E31-3AN

3 leadwire setting for Datascope cable yokes with 2 dummy leads.
SIEMENS/DRÄGER leadwire, 0.95 m (for E202-5XXX cable yokes)



● E32-3IN ● E32-3AN



● E41-5I ● E41-5A

GE leadwire, 0.95 m (for E204-5XXX cable yokes)



● E42-5I ● E42-5A



● E51-3I ● E51-3A

HP/PHILIPS leadwire, 0.95 m (for E205-3XXX cable yokes)



● E52-3I ● E52-3A



● E51-5I ● E51-5A

HP/PHILIPS leadwire, 0.95 m (for E205-5XXX cable yokes)



● E52-5I ● E52-5A



● E61-5I ● E61-5A

SPACELABS leadwire, 0.70 m (E206-5XXX cable yokes)



● E62-5I ● E62-5A



● E61-3I ● E61-3A

SPACELABS leadwire, 0.70 m (for E206-3XXX cable yokes)



● E62-3I ● E62-3A



● E91-5I ● E91-5A

DATASCOPE leadwire, 0.95 m



● E92-5I ● E92-5A

ECG Cables for Patient Monitoring

Leadwires

▼ Yoke Side

▼ Patient Side



- E110GX/5-OR Grouped Leadwire with ESU protection for Philips/HP M1668A, M1663A, M1949A and FMT E105 cable yokes. (C/R/L/ 100cm, F/N 160cm) *Only Grabber option is available.*



Example:
E110GI/5-OR leadwire with
5 grabber end, IEC Color Code
(N, L, F, R, C/C1)



PHILIPS leadwire set, 0.95 m (for E208-3XXX cable yokes)



PHILIPS leadwire set with ESU protection, 0.95 m (for E208-3XXX cable yokes)

Only Grabber.



GE leadwire set, 0.95 m (for E204 cable yokes)



GE leadwire set, 0.95 m (for E204 cable yokes)

● IEC color code ● AHA color code

Single Pin Lead Sets for MultiMed plus (MS20093) and MultiMed plus OR (MS20094) multiparameter cables

For use with Drager Infinity M300 / M540 patient monitors

- **MT32I** 3 lead, IEC, 1 m (OEM P/N MS16159)
- **MT32A** 3 lead, AHA, 1 m (OEM P/N MS16231)



- **MT52I** 5 lead, IEC, 1 m (OEM P/N MS16158)
- **MT52A** 5 lead, AHA, 1 m (OEM P/N MS16546)



- **MT62I** 6 lead, IEC, 1 m (OEM P/N MS16157)
- **MT62A** 6 lead, AHA, 1 m (OEM P/N MS16547)



Multi Parameter Cables

- **FMT-MTM/01** (OEM P/N # 3368391) Multi Parameter Cable is designed to monitor 3 or 5 lead ECG, impedance respiration, temperature and SpO2.

FMT-MTM/01 Multi Parameter Cables are compatible with SC6000, SC7000, SC8000, SC9000, Drager Infinity, Delta and Delta XL monitors (with Nellcor SpO2 option).



Holter ECG Leadwires



- **HL01-7** Holter leadwire, 7 leadwire set, 1.6 mm safety DIN, snap patient end, 90 cm, AHA.

- **HL01-10** Holter leadwire, 10 leadwire set, 1.6 mm safety DIN, snap patient end, 90 cm, AHA.



- **HL02** Mortara holter leadwire, 10 leadwire set, snap patient end, IEC.



- **HL03A** Reynolds Medical ST80/VS20 Event Recorder single channel holter leadwire, 2 lead, snap patient end, IEC.

- **HL03B** Reynolds Medical ST80/VS20 Event Recorder dual channel holter leadwire, 3 lead, snap patient end, IEC.

- **HL10-XX** Biomedical Instruments holter leadwire, 5/7/10 leadwire set, snap patient end, IEC.



- **HL11-XX** DMS holter leadwire, 5/7/10 leadwire set, snap patient end, IEC.

FMT Disposable ECG Leadwires

- FMT Disposable ECG Leadwires reduce the risk of cross contamination
- The Cables are available in 3 & 5 leadwire forms
- FMT disposable ECG leadwires provide compatibility with many of the OEM in the market.




PHILIPS leadwire set (for E208 cable yokes)



GE leadwire set (for E204 cable yokes)

FMT is registered trademark of **Metko Ltd.**



Headquarter & Factory

İvedik O.S.B. Ağaç İşleri Sanayi Sit.
1354. Cad. 1358.Sok. No:9
06378 Yenimahalle

ANKARA / TÜRKİYE

Tel (pbx) : +90 312 387 12 46
Fax : +90 312 387 12 51

Branch

Atatürk Mah. Sedef Cad. 38.Ada
Ata Plaza 3-5 Blokları Kat:2 No:57
34758 Ataşehir

İSTANBUL / TÜRKİYE

Tel (pbx) : +90 216 455 91 21
Fax : +90 216 455 48 24

www.metkold.com

Other products, logos and company names mentioned herein may be trademarks of their respective owners and are not to be taken as an endorsement of affiliation with FMT®. The information in this catalogue has been carefully checked and it is believed to be accurate. Images in this catalogue are representative of items offered. Actual product may have slight visual differences. In the interest of continued product development, Metko reserves the right to make changes and improvements to this catalogue and the product it described any time, without notice or obligation.